

Sidoarjo, 30 Maret 2026

Kepada Yth,
Customer : Rumah Sakit Prima Husada
Address : Jl. Raya Mondoroko, Mondoroko, Banjararum, Kec. Singosari,
Kabupaten Malang, Jawa Timur 65153
Up : Direktur

Dengan hormat,

Kami dari PT. Harlisa Anugerah Sejati selaku sole agent di Indonesia. Yang bergerak dalam bidang critical care, physiotherapy, rehabilitation, orthopedics, ophthalmology dan instruments dengan ini menyampaikan penawaran harga sebagai berikut :

No	Product Name	Brand	Price	Discount	Total Price
1	EEG Type F 24 Channels Standard	NCC	Rp. 250.326.647	Rp. 75.097.994	Rp. 175.228.653
2	EEG Type F 32 Channels with Trolley	NCC	Rp. 475.620.629	Rp. 142.686.189	Rp. 332.934.441
3	EEG Type D 24 Channels Standard	NCC	Rp. 293.325.467	Rp. 87.997.640	Rp. 205.327.827
4	EEG/PSG Type D 36 Channels with Accessories	NCC	Rp. 586.650.934	Rp. 175.995.280	Rp. 410.655.653
5	EEG/PSG/ERP Type D 48 Channels with Accessories and Trolley	NCC	Rp. 1.150.179.816	Rp. 345.053.944	Rp. 805.125.872
6	EEG/PSG/ERP Type D 60 Channels with Accessories and Trolley	NCC	Rp. 1.461.134.354	Rp. 438.340.306	Rp. 1.022.794.048

- Note

➤ Masa berlaku penawaran

➤ Sistem Pembayaran
- : - Harga sudah termasuk instalasi, uji fungsi, dan pelatihan operator

- Barang Ready Stok 10 Hari

- Indent 2 – 3 Bulan

- Harga Sudah Termasuk Ongkos Kirim

- Harga Sudah Termasuk PPN

: 30 hari

: - DP 1 50% : Saat Terbit PO

- Pelunasan 50% : Setelah Instalansi

Demikian penawaran ini kami sampaikan dan kami tunggu kabar baiknya untuk pesanan barang tersebut. Atas perhatian dan kerjasamanya kami ucapkan terima kasih


Hormat kami,
PT. Harlisa Anugerah Sejati



Alice Setya Wati
Manager of Sales Departmen

Routine EEG/PSG/ERP System

Type D EEG



Key features

- 24/36/48/60/72 channels available
- Fiber optical transmission and DC battery power operation used to prevent interference from AC
- Built-in impedance testing function to check electrodes wear situation
- High accuracy EEG with 24 bit precision
- Support EDF/TXT formats
- Ergonomically designed single shielded cup/clip electrodes with DIN interface
- Expandable PSG recording: EOG, Air Flow, Snoring, ECG, Spo2, thoracic and abdominal respiration, synchro-nized EEG examination together with PSG enabling for more sophisticated clinical applications
- Optional Acoustic & Visual stimulator for ERP examination to check patient's recognition situation, including P300, MMN, CNV and P50 etc
- Optional HD infra-red camera enabling simultaneous acquisition, editing, and rapid display of video signals and EEG signals

EEG System configuration

» Various EEG Electrode



General bracket electrode Cap



Bracket silver electrode



Silver ear electrode



Single bracket electrode cable




Integrated bracket electrode cable



Single cup electrode cable



Integrated cup electrode cable



Integrated cap electrode cable

» PSG Electrode



EOG/EMG



Air Flow



Spo2



Snoring



ECG




Respiration

» Optional configuration



Photoc Stimulator



Infra-red Camera



ERP Accessories

Clinical application

- Discrimination between cerebral organic pathologies and functional disorders
- Diagnosis of diverse cerebral ailments
- Localization assessment of intracranial lesions
- Evaluation of therapeutic effectiveness and prognostication for diseases
- Assessment of disease progression and cerebral functional status
- Utilization cerebral injuries and epileptic manifestations

For more, please contact us:

Diagnostic: EEG, EMG

Surgical monitoring: Cynapse IONM, Smart IONM

Rehabilitation: Biofeedback system, Micro-current stimulator, Surface EMG system



Distributed by
PT Harlisa Anugerah Sejati
JL. Sriti Blok HH No.15, RT 002/RW 003, Pabean
Sedati, Sidoarjo, 61253, Jawa Timur - Indonesia
Tel: +62 31 9968 3467, www.harlisa.co.id
E-mail: info@harlisa.co.id

Routine/Ambulatory EEG System

Type F EEG

Key features

- 24/32 Channels available
- Compact design, portable and easy to carry
- Wi-Fi transmission facilitating real-time recording of EEG waveforms, allowing patients to move freely
- Patient's marker button utilized to indicate the occurrence of epileptic seizures during routine or long-term EEG recordings
- Low power consumption supporting up to 24 hours EEG waveform recording with external rechargeable batteries
- Dynamic recording with 16GB flash memory card, capable of supporting 4-5 days case storage and enabling long term EEG
- All-in-one electrode cable for easy connection
- DC Power supply eliminating risk of electric shock to patients and preventing signal interference
- Enabling 24/32 CH EEG or 16/24CH+8 Bipolar channels for ECG and EMG
- Optional HD infra-red camera enabling simultaneous acquisition, editing, and rapid display of video signals and EEG signals



Type F EEG-Holter



5
Biocompatibility test



6
Accelerated aging test



8
Simulated transportation test



14
Electromagnetic compatibility (EMC) test



21
Performance indicator test



385
Electrical safety monitoring



880
Software test and detection



Type F EEG-Routine

Key features

- 16/24/32/64 Channels available
- Compact design, portable and easy to carry
- Wi-Fi transmission facilitating real-time recording of EEG waveforms, allowing patients to move freely
- The low-power consumption amplifier capable of supporting remote recording for up to 72 hours with dry batteries
- All-in-one electrode cable for easy connection
- Built-in impedance testing function to check electrodes wear situation

Clinical application

- Discrimination between cerebral organic pathologies and functional disorders
- Diagnosis of diverse cerebral ailments
- Localization assessment of intracranial lesions
- Evaluation of therapeutic effectiveness and prognostication for diseases
- Assessment of disease progression and cerebral functional status
- Screening for neonatal cerebral injuries and epileptic manifestations
- Utilization cerebral injuries and epileptic manifestations

For more, please contact us:

Diagnostic: EEG, EMG

Surgical monitoring: Cynapse IONM, Smart IONM

Rehabilitation: Biofeedback system, Micro-current stimulator, Surface EMG system



Distributed by
PT Harlisa Anugerah Sejati
 JL. Sriti Blok HH No.15, RT 002/RW 003, Pabean
 Sedati, Sidoarjo, 61253, Jawa Timur - Indonesia
 Tel: +62 31 9968 3467, www.harlisa.co.id
 E-mail: info@harlisa.co.id



PT. PETAN DAYA MEDICA

APL TOWER
Podomoro City
Lantai 12, Unit 2
Jl. Letjen S. Parman Kav. 12
JAKARTA 11470
Tlp: 29339259; 29339260
Fax: 29339258
Email : petandaya@yahoo.com.sg

No. 2965/PNW/PDM-JKT/XII/2025

Jakarta, 11 Desember 2025

Hal : **Penawaran Peralatan EEG 32 Channel**

Kepada Yth.

Direktur Utama

RS Prima Husada Sukorejo

Jl. Raya Surabaya- Malang

Pasuruan

Jawa Timur

Dengan Hormat,

Bersama ini perkenalkan kami selaku agen tunggal di Indonesia dari **"Natus" (USA)**, untuk menawarkan peralatan neuro diagnostik, dengan rincian sebagai berikut:

- | | |
|---------------------|---------------------------|
| 1. Nama Alat | : EEG 32 Channel |
| Merek | : Natus (USA) |
| Type | : Nicolet |
| Harga | : Rp 795.000.000,- |

Kondisi Penawaran

- | | |
|---------------------------------|--|
| Harga | : Belum Termasuk PPN 11% |
| Spesifikasi & Brosur | : Terlampir |
| Pembayaran | : Uang muka 40% saat P.O kami terima, 60% sisanya pada saat barang diserahkan |
| Penyerahan Barang | : Ready Stock (selama persediaan masih ada)/Indent 60–90 hari |
| Garansi | : 1 (satu) tahun berikut suku cadang dan Jasa pemeliharaan |
| Training | : Operator / perawat dilakukan di tempat |
| Penawaran | : Berlaku 30 hari dari tanggal Penawaran |
| | Harga belum termasuk Ongkos Kirim |
| | Tidak berlaku jika USD 1,- melebihi Rp 16.700,- |

Demikian penawaran ini kami sampaikan, atas kepercayaan dan kesempatan yang sudah diberikan kami ucapkan terima kasih.

Hormat kami,

PT. Petan Daya Medica

A/N

Andrew Djuli

Manager



PT. PETAN DAYA MEDICA

APL Tower
Lantai 12 Suite 2
Padomoro City
Jl. Letjen S. Parman Kav. 28
Jakarta 11470
Phone :
(021) 293 392 59 / 293 392 60
Fax :
(021) 293 392 58
Email :
petandaya@yahoo.com.sg

Nama Alat : Electro Encephalograph / EEG
Merek : Natus
Type : Nicolet EEG
Negara Asal : USA

Description	Specification	QTY
		Set
Operating Environment (in use)	: Temperature: 15.6 to 32.2° C, (60 to 90° F)	
	Relative Humidity: 20-80%, non-condensing	
	Altitude: 0-3km, (0-10,000 ft)	
Non-Operating Environment (in storage)	: Temperature: 17.7 to 55° C, (O to 132° F)	
	Relative Humidity: 10-90%, non-condensing	
	Altitude: 0-12km, (0-40,000 ft)	
Desktop Computer (with local trolley)	: CPU Intel i5 or better with 8GB RAM (min.)	
	Operating System Microsoft® Windows 10	
	User Interface Keyboard and mouse,	
	Storage Devices: SSD 256 GB with 1 TB and HDD External 2TB	
	DVD+R/W Drive	
Monitor/Display	: Pixel Resolution 1280 X 1024	
	Min 24 " min. wide aspect LCD (min) can be upgraded upon necessary	
VEEG software (with v32)	: EEG software	
	Upgradeable with Brain Mapping License Support	
	EEG aquisition &review software	
	Feature:	
	Integrated Impedance display for see your inpedance value right at the bedside	
	Integrated SPO2 for added value with SPO2 enabling sleep monitoring	
	Etherned amplifier interface for industry standard conection	
	9 auxiliary pair(AC or DC) for ability to connect devices such as respiratory belts, thermisters, ect	
	One High level DC input(v32 only) for ability to connect various devices such as CPAP	
	12 High level DC inputs (v44 only for ability to connect multiple devices such as CPAP	
Report generation software	: Topographical Brain Mapping (Optional)	
	o Amplitude Mapping	
	o Frequency Band Power mapping	
	o Coherence Mapping	
EEG Display	: Sec/Page 2, 5, 10, 20, 30, 60, 120, 240, 300, 600, 1200	
	mm/Sec 6, 8, 10, 15, 30, 60, 120, 240	



PT. PETAN DAYA MEDICA

APL Tower
Lantai 12 Suite 2
Podomoro City
Jl. Letjen S. Parman Kav. 28
Jakarta 11470
Phone :
(021) 293 392 59 / 293 392 60
Fax :
(021) 293 392 58
Email :
petandaya@yahoo.com.sg

Nama Alat : Electro Encephalograph / EEG
Merek : Natus
Type : Nicolet EEG
Negara Asal : USA

Description	Specification	QTY
		Set
Sensitivity	: 10, 20, 30, 50, 70, 100, 150, 200, 300, 500, 700, 1000, 2000, 5000 μ V/cm	
	1, 2, 3, 5, 7, 10, 15, 20, 30, 50, 70, 100, 200, 500 V/m μ V	
	High Filters 15, 30, 35, 40, 50, 60, 70, 100, 150, 200, 300, 500, 1000, 1500 Hz	
	Low Filters Off 0.016, 0.16 0.3, 0.5, 1, 1.6, 2, 3, 5 Hz	
	0.03,0.05, 0.1,0.15,0.3, 0.5,1,10 seconds	
	Notch Filter Off, 50/60 Hz	
	Photic stimulator LED or flash	
Print out	: Printer	
	Network	
	10/100/1000 Mb ethernet (standard)	
Amplifier	: System Configurations	
	ADC Resolution Voltage = 0.153 μ V	
	DC Offset Tolerance + 900 mV	
	Channels (Inputs) 32 EEG (9 off configurable as bipolar AC) ,1 off configurable as DC	
	Maximum Input Range $\pm$ 5 mV	
	Bandwidth 0.053 - 500 Hz	
	Noise < 1.5 μ V pk-pk @ 0.1 - 100 Hz	
	Input Impedance > 100 MW (common mode)	
	CMRR at Patient Inputs > 115 dB @ 50 - 60 Hz, with active patient ground connected	
	Channel Crosstalk < -40 dB	
	Amplifier Sample Rate (under software control)	
	125, 250, 500, 1000, 2000	
	Calibration Square wave, 1, 5, 10, 20 sec period,	
	10, 50, 100, 1000 μ V amplitude	
	Input Bias Current < 5 nA	
	Anti-Aliasing Filter Cut Off Frequency 500 Hz	
	Differential Input Impedance 40 M Ω	
	Interface to Amplifier Ethernet	
	Channel Hardware Gain 410	
	Deblock Yes	
	Auxiliary inputs	
	1Hi-level, non isolated inputs for connection of external devices(e.g	
	Spo2, Co2 monitor, etc)	
	•Analog/Digital Converter 16 bits	
	•Maximum Input Range $\pm$ 2.5V	
	•ADC Resolution 76.3 μ V	
	•Bandwidth DC - 500 Hz	

**POPULASI ELECTROENCEPHALOGRAPH**

No	Type EEG	UNIT	CUSTOMER
1	Nicolet vEEG 32 Ch	1	RS PGI Cikini Jakarta
2	Nicolet vEEG 32 Ch	1	Klinik Zam Zam (dr. Nizmah) Jakarta
3	Nicolet vEEG 32 Ch	1	Lab Pramita Matraman Jakarta
4	Nicolet vEEG 32 Ch	1	RS. Atmajaya Jakarta
5	Nicolet vEEG32 Ch	1	RS. MMC Jakarta
6	Nicolet vEEG 44 Ch	1	RS Premier Jatinegara Jakarta
7	Nicolet vEEG 44 Ch	1	RSUP.Pertamina Jakarta
8	Nicolet vEEG 32 Ch	1	RS. Fatmawati Jakarta
9	Nicolet vEEG 32 Ch	1	RSUD Tarakan Jakarta
10	Nicolet vEEG 32 Ch	1	RS. Puri Cinere Jakarta
11	Nicolet vEEG 32 Ch	1	RS. Permata Cibubur Jakarta
12	Nicolet vEEG 32 Ch	1	RS. Adiyaksa Jakarta
13	Nicolet vEEG 32 Ch	1	RS. Siloam Kebon Jeruk Jakarta
14	Nicolet vEEG 32 Ch	1	RS. Siloam TB. Simatupang Jakarta
15	Nicolet vEEG 32 Ch	1	RSU. Bunda Menteng Jakarta
16	Nicolet wEEG 64 Ch	1	RS. POLRI Jakarta
17	Nicolet vEEG 32 Ch	1	Balai Kesehatan Penerbangan
18	Nicolet vEEG 32 Ch	1	Balai Kesehatan Penerbangan
19	Nicolet vEEG 32 Ch	1	RS. Hermina Depok
20	Nicolet vEEG 32 Ch	1	Lab Parahita Bandung
21	Nicolet vEEG 32 Ch	1	Klinik Samya (dr. Suryani) Bandung
22	Nicolet vEEG 32 Ch	1	Lab Pramita Cirebon
23	Nicolet vEEG 32 Ch	1	Lab Pramita Riau Bandung
24	Nicolet vEEG 32 Ch	1	RS. Majalaya
25	Nicolet vEEG 32 Ch	1	RS. Boromeus Bandung
26	Nicolet vEEG 32 Ch	1	RSUD. Tangerang Selatan
27	Nicolet vEEG 32 Ch	1	RSUD Dr. Drajat Prawiranegara Serang
28	Nicolet vEEG 32 Ch	1	RSU Banten
29	Nicolet vEEG 32Ch	1	RSUD Badung Bali
30	Nicolet vEEG 32Ch	1	RS. Prima Medica-Bali
31	Nicolet VEEG 32Ch	1	RS. Sanglah (Wing) Bali
32	Nicolet VEEG 32Ch	1	RS. Sanglah (Univ) Bali



No	Type EEG	UNIT	CUSTOMER
33	Nicolet VEEG 32Ch	1	RS. Sanglah (Sleep) Bali
34	Nicolet VEEG 32Ch	1	RSUD. Bangli Bali
35	Nicolet VEEG 32Ch	1	RS. Siloam Bali
36	Nicolet VEEG 32Ch	1	RSUD. Sanjiwani Gianyar Bali
37	Nicolet VEEG 32Ch	1	Klungkung Bali
38	Nicolet VEEG 32 Ch	1	RS. Adam Malik Medan SUMUT
39	NicoletOne VEEG 44 Ch	1	RSU. Kaban Jahe SUMUT
40	NicoletOne VEEG 32 Ch	1	RS. Sipirok SUMUT
41	NicoletOne VEEG 32 Ch	1	RS. Bunda Thamrin Medan SUMUT
42	Nicolet vEEG 32 Ch	1	Lab Pramita Medan SUMUT
43	Nicolet vEEG 32 Ch	1	RSUD Binjai SUMUT
44	Nicolet wEEG32 Ch	1	RS. Haji Medan SUMUT
45	Nicolet vEEG 32 Ch	1	RS. Haji Medan SUMUT
46	Nicolet WEEG 32 Ch	1	Klinik GMC Medan SUMUT
47	Nicolet vEEG 32 Ch	1	RS. Siloam Medan SUMUT
48	Nicolet VEEG 32 Ch	1	RS. Royal Prima Medan SUMUT
49	Nicolet vEEG 32 Ch	1	RS.USU Medan SUMUT
50	Nicolet vEEG 32 Ch	1	RSUD. Padang Sidempuan SUMUT
51	Nicolet VEEG 44 Ch	1	RSUD. Sultan Sulaiman Serdang Bedagai SUMUT
52	Nicolet BFM	1	RSUD Tanjung Balai SUMUT
53	Nicolet vEEG 32 Ch	1	RSUD Kisaran SUMUT
54	Nicolet vEEG 32 Ch	1	RSU. Dr. M Djamil Padang
55	Nicolet vEEG 32 Ch	1	Lab Pramita Padang
56	Nicolet vEEG 32 Ch	1	RSUD. Dr. Adnan WD Payakumbuh
57	Nicolet vEEG 32 Ch	1	RS. Siloam Jambi
58	Nicolet vEEG 32 Ch	1	RSUD. Ahmad Ripin Jambi
59	Nicolet vEEG 32 Ch	1	RSUD Rejang Lebong Bengkulu
60	Nicolet vEEG 32 Ch	1	RSJ Lampung
61	Nicolet vEEG 32 Ch	1	RS. Prima Medika Pekanbaru
62	Nicolet EEG v44 Ch	1	National Hospital Surabaya
63	Nicolet vEEG 32 Ch	1	Lab Pramita Surabaya Adityawarman
64	Nicolet VEEG 32 Ch	1	Lab Pramita Surabaya Jemur Andayani
65	Nicolet vEEG 32 Ch	1	Lab Sima Malang
66	Nicolet v44 Ch LTM	1	Lab Pramita Madiun
67	Nicolet vEEG 32 Ch	1	RSUD. Dr. Haryoto Lumajang
68	Nicolet vEEG 32 Ch	1	RSUD Soetomo Surabaya



No	Type EEG	UNIT	CUSTOMER
69	Nicolet wEEG 64 Ch	1	RSUD Soetomo Surabaya
70	Nicolet vEEG 32 Ch	1	RS. Haji Surabaya
71	Nicolet vEEG 32 Ch	1	RSUD. Ngawi
72	Nicolet vEEG 32 Ch	1	RSUD Ponorogo
73	Nicolet vEEG 32 Ch	1	RS. Zainal abidin Banda Aceh
74	Nicolet BFM	1	RS. Zainal abidin Banda Aceh
75	Nicolet vEEG 32 Ch	1	RS. Datu Beru Takengon Aceh Tengah
76	Nicolet vEEG 32 Ch	1	RS. Siloam Yogyakarta
77	Nicolet vEEG 32 Ch	1	Pramita Yogya
78	Nicolet vEEG 32 Ch	1	Lab Pramita Semarang
79	Nicolet wEEG32 Ch	1	RS. Jiwa Solo
80	Nicolet vEEG 32 Ch	1	RS. Sitti Khodijah Palembang
81	Nicolet vEEG 32 Ch	1	RS. M.Husain, Palembang
82	Nicolet BFM	1	RS. M.Husain, Palembang
83	Nicolet vEEG 32 Ch	1	RS. Siloam Palembang
84	Nicolet vEEG 32 Ch	1	RS. AR Bunda Lubuk Linggau
85	Nicolet vEEG 32 Ch	1	RS. YK Madina Palembang
86	Nicolet vEEG 32 Ch	1	RS. Siloam Makassar
87	Nicolet vEEG 32 Ch	1	Klinik dr. Amirudin Aliah Makassar
88	Nicolet vEEG 32 Ch	1	Klinik dr. Herlyani (Manado)
89	Nicolet vEEG 32 Ch	1	RSUD. Kandau Manado
90	Nicolet vEEG 32 Ch	1	RSUD Lattemallama Sopeng
91	Nicolet vEEG 32 Ch	1	RSUD. Tani dan Nelayan Boalemo Gorontalo
92	Nicolet vEEG 32 Ch	1	RSUD Aloe Saboe Gorontalo
93	Nicolet vEEG 32 Ch	1	RSUD Soedarso Pontianak Kalimantan Barat
94	Nicolet vEEG 32 Ch	1	RSUD. Sintang Kalimantan Barat
95	Nicolet vEEG 32 Ch	1	RSJ Singkawang Kalimantan Barat
96	Nicolet vEEG 32 Ch	1	RS. Mitra Medika Pontianak Kalimantan Barat
97	Nicolet vEEG 32 Ch	1	RSUD Tarakan kalimantan Utara
98	Nicolet vEEG 32 Ch	1	RSUD Pangkalan Bun Kalimantan Tengah
99	Nicolet vEEG 32 Ch	1	RS. Siloam Balikpapan Kalimantan Timur
100	Nicolet vEEG 32 Ch	1	RS. Pertamina Balikpapan Kalimantan Timur
101	Nicolet vEEG 32 Ch	1	RS. Bhayangkara Samarinda Kalimantan Timur
102	Nicolet vEEG 32 Ch	1	RSUD. Kota Mataram NTB
103	Nicolet vEEG 32 Ch	1	RS. Siloam Kupang NTT
104	Nicolet vEEG 32 Ch	1	RS PON Jakarta
105	Nicolet vEEG 32 Ch	1	RSUD Dr. Saiful Anwar, Malang Jawa Timur
106	Nicolet vEEG 32 Ch	1	RSUD Belitung Timur, Bangka Belitung
107	Nicolet vEEG 32 Ch	1	RSUD dr. Chasbullah Abdul Madjid, Bekasi, Jawa Barat



No	Type EEG	UNIT	CUSTOMER
108	Nicolet vEEG 32 Ch	1	RSUD Bung Karno, Surakarta, Jawa Tengah
109	Nicolet vEEG 32 Ch	1	RSUD Gambiran, Kediri, Jawa Timur
110	Nicolet vEEG 32 Ch	1	RSUD Kabupaten Pidie Jaya, Aceh
111	Nicolet vEEG 32 Ch	1	RS dr. Abdul Radjak, DKI Jakarta
112	Nicolet vEEG 32 Ch	1	RSUD Tugu Koja, Jl. Walang Permai , DKI Jakarta
113	NicoletOne VEEG 44 Ch	1	RS Premier Bintaro, Jl. MH Thamrin No. 1, Tangerang
114	Nicolet vEEG 32 Ch	1	RSUD Sapiro, Jl. Rumah Sakit No. 1, Tapanuli, Sumatera Utara
115	EEG 64 Ch	1	RSUD Kota Malang, Jl. Rajasa No.27, Malang, Jawa Timur
116	EEG 64 Ch	1	RSUD Dr. Iskak Tulungagung ,Jl. Dr. Wahidin Sudiro Husodo, Tulungagung, Jawa Timur
117	Nicolet vEEG 32 Ch	1	RSUD Padangan, Jl. Dr. Soetomo No. 02, Padangan Jawa Timur
118	EEG 64 Ch	1	RSUD Koja, Jl. Deli No.4, Jakarta Utara
119	EEG 64 Ch	1	Balai Kesehatan Penerbangan, Jl. Kota Baru Bandar Kemayoran Blok B11 Kav. No. 4, Jakarta
120	EEG 64 Ch	1	RSU Dr. Cipto Mangunkusumo, Jl. Diponegoro No. 71, Jakarta Pusat
121	EEG 64 Ch	1	RS Grha Kedoya, Jl. Panjang No.26, Jakarta Barat
122	Nicolet EEG 32 Ch	1	Siloam Hospitals Manado,Jl. Sam Ratulangi No.22, Manado, Sulawesi Utara
123	Nicolet EEG 32 Ch	1	RS Bhayangkara Makassar, Jl. Andi Mappaodang No.63, Makassar, Sulawesi Selatan
124	Nicolet EEG 32 Ch	1	RSU St. Madyang Palopo, Jl. Andi Kambo No.87, Palopo, Sulawesi Selatan
125	EEG 64 Ch	1	RSUD Kota Cilegon, Jl. Kapten Pierre Tendean km.03, Cilegon, Banten
126	EEG 64 Ch	1	RS Mitra Keluarga Satelit, Jl. Satelit Indah II, Surabaya, Jawa Timur
127	Nicolet EEG 32 Ch	1	RS Awal Bros Panam Pekanbaru, Jl. HR. Soebrantas No.88, Pekanbaru, Riau
128	Nicolet EEG 32 Ch	1	Sentra Medika Hospital Minahasa Utara, Jl.Ir. Soekarno, Minahasa Utara, Sulawesi Utara
129	Nicolet EEG 32 Ch	1	RSUD Dr H Abdul Moeloek, Jl. Dr. Rivai No.6 Bandar Lampung, Lampung
130	EEG 64 Ch	1	RSUD Dr H Abdul Moeloek, Jl. Dr. Rivai No.6 Bandar Lampung, Lampung
131	Nicolet EEG 32 Ch	1	Siloam Hospitals Purwakarta , Jl. Raya Bungursari No.1, Purwakarta, Jawa Barat
132	EEG 64 Ch	1	RSUD Aceh Tamiang, Jl. Kesehatan No.1 Kabupaten Aceh Tamiang, Aceh
133	EEG 64 Ch	1	RSUD Soedarso Pontianak, Jl. DR. Soedarso No.1 Kalimantan Barat
134	EEG 64 Ch	1	RSUD Munyang Kute Redelong, Jln. Simpang Teritit-Pondok Baru, Aceh



No	Type EEG	UNIT	CUSTOMER
135	EEG 64 Ch	1	RSUD Rantauprapat, Jl. KH. Dewantara No.129 Labuhan Batu, Sumatera Utara
136	EEG 64 Ch	1	RSUD Datu Beru, Jl. Qurata Aini No. 153 Kebayakan, Takengon, Aceh Tengah
137	Nicolet EEG 32 Ch	1	RSUD Bima Jl. Langsung No. 1, Bima, Nusa Tenggara Barat
138	Nicolet EEG 32 Ch	1	RS Al Irsyad Surabaya, Jl. KH. Mas Manysur No. 210-214 Surabaya, Jawa Timur
139	EEG 64 Ch	1	RSUD Kalideres, Jl. Raya Daan Mogot, Jakarta
140	EEG 64 Ch	1	RSUP Dr. Hasan Sadikin, Jl. Pasteur No. 38, Bandung, Jawa Barat
141	Nicolet EEG 32 Ch	1	RSUD Sukamara, Jl. Tjilik Riwut No.km. 5.5, Sukamara, Kalimantan Tengah
142	Nicolet EEG 32 Ch	1	RSUD Raden Mattahter, Jl. Letjen Suprpto No.31, Jambi
143	Nicolet EEG 32 Ch	1	RS Awal Bros Dumai, Jl. Jend. Sudirman No.88, Dumai, Riau
144	EEG 64 Ch	1	RSUPT Vertikal Surabaya, Jl. Indrapura No.17, Surabaya, Jawa Timur
145	EEG 64 Ch	1	RSUD Nurdin Hamzah, Jl. WR. Supratman Kel. Parit Culum I, Jambi
146	EEG 64 Ch	1	RSUD Mas Amsyar, Jl. Tjilik Riwut No.1 Kasongan, Kalimantan Tengah
147	Nicolet EEG 32 Ch	1	RS Awal Bros Batu Aji, Jl. Trembesi Point, Kibing, Batam, Kepulauan Riau
148	Nicolet EEG 32 Ch	1	RS Awal Bros Hangtuah, Jalan Hangtuah, Pekanbaru, Riau
149	EEG 64 Ch	1	RSUD Pasar Rebo, Jl. TB Simatupang No.30, Jakarta Timur
150	Nicolet EEG 32 Ch	1	RSUD Jatisari, Jl. Raya Jatisari No.3, Karawang, Jawa Barat
151	Nicolet EEG 32 Ch	1	RSUD Mamuju, Jl. Kurungan Bassi No.1, Mamuju, Sulawesi Barat
152	Nicolet EEG 32 Ch	1	RSUD Kabupaten Kediri, Jl. Pahlawan Kusuma Bangsa No.01, Kediri, Jawa Timur
153	EEG 64 Ch	1	RSUD Srengat, Jl. Raya Dandong, Blitar, Jawa Timur
154	Nicolet EEG 32 Ch	1	Eka Hospital MT Haryono, Jakarta
155	Nicolet EEG 32 Ch	1	Siloam Hospitals Palangkaraya, Jl. RTA Milono, Palangka Raya
156	Nicolet EEG 32 Ch	1	Siloam Hospitals Lippo Village Karawaci, Jl. Siloam No. 6, Tangerang, Banten
157	Nicolet EEG 32 Ch	1	RSIA Murni Teguh Rosiva, Jl. Bangka no. 15, Sumatera Utara
158	Nicolet EEG 32 Ch	1	RS Prima Medika, Jl. Raya Sesetan No.10, Denpasar, Bali
159	Nicolet EEG 32 Ch	1	RS Awal Bros Ujung Batu, Jl. Jendral Sudirman No.314, Riau

No : 250105/PNWR/PANAMED/XII/2025
Hal : Penawaran Harga EEG Nihon Kohden

Jakarta, 12 Desember 2025

Kepada Yth:
Direktur Utama
Rumah Sakit Prima Husada
Jl. Raya Mondoroko, Mondoroko, Banjararum, Kec. Singosari, Kabupaten Malang
Jawa Timur 65153

Dengan Hormat,

Bersama ini perkenalkan kami PT. Panamed Anugerah Abadi selaku Distributor Alat Kesehatan Nihon Kohden di Indonesia, untuk menawarkan peralatan kami sebagai berikut:


Cat. No	Qty	Nama Barang	Status Barang	Harga (Rp)	Total (Rp)
EEG-1200K	1	Nihon Kohden EEG-1200K 32 +14 Channel with 3D Brain Mapping	Ready Stock	840,000,000	840,000,000
	1	Nihon Kohden HD Video Camera & EEG Video Software		150,000,000	150,000,000

Kondisi Penawaran sebagai berikut:

Spesifikasi & Brosur : Terlampir
Harga : Belum termasuk ongkos kirim dan PPN 11%
Pembayaran : Uang muka 60% saat PO kami terima, 40% pada saat barang diserahkan
Garansi : 1 (satu) tahun garansi pabrik
Training : 1 (satu) operator
Penawaran : Berlaku 45 hari sejak tanggal penawaran
: Ready selama persediaan masih ada

Demikian penawaran ini kami sampaikan, atas kepercayaan dan kesempatan yang telah diberikan kami ucapkan terima kasih.

Hormat kami,



PT. Panamed Anugerah Abadi
Hospital and Laboratory Specialist

Vania
Direktur

Neurofax



Electroencephalograph

EEG-1200J/K



Fighting Disease with Electronics

 **NIHON KOHDEN**



EEG-1200

Leading healthcare



User-friendly

Efficient


Expandable

Accurate



the future
to the future

into the future...




The possibilities are endless.

User-friendly


Customizable new main menu

You can register up to 10 examination protocol buttons per page on the main menu. Each button has user-defined settings for an examination. The settings can also be changed for different examination conditions.



3D Voltage mapping—Fast review and advanced EEG analysis

Whole head maps provide a complete overview and a better interpretation of the topography of EEG abnormalities. Just click on a detected pattern to obtain 3D maps of the whole head. Click on a particular view to obtain a series of maps showing the change over time.




Note window—Simple copying of waveform parts

You can save up to 1,000 sections of waveforms for comparison by dragging and dropping. Up to 100 copied waveforms can be registered as sample data for comparison with other patients.

Zoom window

You can easily magnify the waveforms by dragging them. Amplitude and latency of magnified waveforms can be measured and printed.



Screen comment tags

Up to 100 tags can be attached to an EEG file for later reference or messages to the reviewer. The tags can include Word documents, Excel spreadsheets, images and other files.


DSA trendgraph

Frequency components of EEG and the amplitudes of each frequency are displayed as a DSA (Density Spectral Array) trend graph on the review screen and EEG Scope. The DSA lets you find epileptic seizures of a specific EEG frequency band over a long period of time at a glance.

Efficient


EEG Scope— Data review during acquisition

Comparison Mode: The EEG Scope function lets you look back into and review previous epochs of EEG while simultaneously monitoring current EEG acquisition. One side of a split screen shows previous EEG epochs and the other side shows the current EEG.



Advanced EEG report generation— Using NeuroReport (standard)

NeuroReport has various item templates for reports and you can create customized reports for different examination types, physicians or other criteria. You can also create an EMU examination report which has screen comments automatically assigned to seizure events. Reports are saved in a database and you can quickly search for reports by patient information, measurement values and other criteria.



Junction box: JE-921A

All-in-one solution: One junction box for routine EEG and PSG.

This advanced junction box integrates 32 channel EEG input and SpO₂/ETCO₂ inputs. JE-921A provides the highest signal quality and maximum reliability. The QI-122A input box converter connects JE-921A to the EEG over a LAN network. This provides LAN connectivity and expandability for an EEG system in the EMU and Sleep center.



- Extra PSG sensor capability—
Input jacks for analog signals from external instruments such as CPAP
- Unique technology—Built-in SpO₂/ETCO₂ inputs
- 4 channels DC input



Expandable

256 channels, JE-120A junction box

JE-120A junction box provides up to 256 channels. It can be used to determine focus of epilepsy by using grid electrode and depth electrodes. Maximum 10,000 Hz high sampling rate enables to measure High Frequency Oscillation (HFO) and up to 10 seconds of Time Constant makes DC shift visible.



- 24 bit/10 kHz high sampling rate
- SpO₂/ETCO₂ measurement is available
- DC 16 inputs
- Up to 4 pairs of bipolar (when using JE-125AK)

JE-125AK or JE-225AK mini flat junction box
(1 to 64 ch)

JE-226AK mini flat junction box
(65 to 128 ch)

JE-227AK mini flat junction box
(129 to 192 ch)

JE-228AK mini flat junction box
(193 to 256 ch)


EEG input over a LAN network

QI-123A input box converter connects junction box to the EEG over a LAN network. This provides LAN connectivity and expandability for an EEG system in the EMU.



New Sophisticated system for functional Brain Mapping Test

Sophisticated system of switch box PE-210AK and extension unit MS-120BK with JE-120A amplifier enhances the efficiency of brain functional mapping tests, which are conducted during intracranial EEG monitoring for a patient with epilepsy. This system enables to select by software control active and reference electrodes for stimulation and to generate a complete stimulation session report.




Accurate

EEG Trendprogram software: A variety of EEG trendgraph QP-160A

Synchronized digital video for EEG systems


- Amplitude-integrated EEG
- FFT analysis such as DSA, FFT Power Ratio, etc
- Burst Suppression Ratio
- DC trendgraphs such as SpO₂, ETCO₂, and HR



EEG Autoreport software, QP-270A

Automatic report of routine EEG with interpretation of main background activity

- Artifacts detection: blink, horizontal eye movement, EMG, earlobe activation, electrode pop
- EEG pattern analysis: Background activity, Spike, Photic activation
- Designed by Prof. Hiroshi Shibasaki and Prof. Masatoshi Nakamura



Digital video software: Video-link QP-110AK

Synchronized digital video for EEG systems

- Patient images synchronized with the EEG waveforms can be recorded
- Video clip and snapshot functions are available
- IP Camera for High-Definition video is now available



Spike detector software: Spike Detector QP-251AK

On-line and off-line spike and seizure detection vastly improved with greater accuracy

Analysis software: EEGFocus QP-211AK


Review, mapping, remontaging, filtering and FFT

- Voltage and current source density mapping
- Automatic rejection of eye blink artifact without modifying the original data
- Automatic detection and averaging of similar waves
- Source imaging

EEG mapping software: QP-220AK


Real-time and basic EEG mapping software

- Real-time and off-line mapping
- Up to eight frequency maps (seven power/voltage maps at seven different frequency bands and one map of all frequency bands)
- Power/Voltage spectra for up to 32 channels of EEG waveform data
- Edge frequency, average frequency, median frequency or peak frequency for each spectrum is indicated with a mark



Sleep analysis software: POLYSMITH Nihon Kohden Station, PS-ONLINE

- Display, store and analyze sleep data
- Reformatting of individual channel filters and sensitivities or montage
- Customize the system to best suit your needs
- Time-link trendgraphs superimposed on waveform data with epoch details



Composition examples



Routine EEG composition

Electroencephalograph	EEG-1200J/K
PC unit	CC-120AJ/AK
Electrode junction box	JE-921A
Flash lamp assembly	LS-703A
Photoc stimulator control unit	LS-120AJ/AK
Cart	KE-122A
Stand	KC-001A
LCD display	local purchase



EMU Composition

Electroencephalograph	EEG-1200J/K
PC unit	CC-120AJ/AK
256 channel electrode junction box	JE-120A-256
Input box converter	QI-123A
Mini flat junction box	JE-125AK
Mini flat junction box	JE-226AK
Mini flat junction box	JE-227AK
Mini flat junction box	JE-228AK
Flash lamp assembly	LS-703A
Photoc stimulator control unit	LS-120AJ/AK
Cart with printer table	KD-029A
Stand	KC-001A
LCD display	Local purchase
Extension unit	MS-120BK
Switch box	PE-210AK

Model Suffixes

EEG-1200 has the following suffixes:

J: 110-127 V AC operation

K: 220-240 V AC operation

Spike Detector software is a product of Persyst Development Corporation.

EEGFocus is a product of BESA GmbH.

Polysmith is a trademark of Neurotronics Inc.

Windows and the Windows logo are registered trademarks of Microsoft Corporation.

This brochure may be revised or replaced by Nihon Kohden at any time without notice.



NIHON KOHDEN CORPORATION

1-31-4 Nishiiochiai, Shinjuku-ku, Tokyo 161-8560, Japan
Phone +81 (3) 5996-8036 Fax +81 (3) 5996-8100

www.nihonkohden.com