

Surabaya, 27 Maret 2026

No : EI/074.03/RSPH/CTR/2026
Perihal : Penawaran Endoscopy
Lampiran : 1 (satu) bundle

Kepada Yth.

Direktur

RS Prima Husada Malang

Jl. Raya Mondoroko, Mondoroko, Banjararum

Kec. Singosari, Kabupaten Malang, Jawa Timur 65153

Dengan hormat,

Kami dari PT. Endo Indonesia selaku Distributor/Sole Agent dari peralatan kesehatan, Maka dalam kesempatan ini, kami ingin menyampaikan penawaran **Endoscopy** sebagai berikut:

NO	DESCRIPTION	QTY	MERK	HARGA PENAWARAN	HARGA SETELAH DISC 30%
	ENDOSCOPY SYSTEM				
1	Video Endoscope System EI.EC.v3	1 set	ENDO - Sonoscape	Rp 4,398,771,429	Rp 3,079,140,000
2	Video Endoscope System EI.EC.v2	1 set	ENDO - Sonoscape	Rp 3,426,757,143	Rp 2,398,730,000
3	Video Endoscope System EI.EC.v7	1 set	ENDO - Sonoscape	Rp 2,214,285,714	Rp 1,550,000,000

- Harga sudah termasuk:

1. PPN 11%

- Ketentuan Sistem Pembayaran :

1. Stok : Ready

2. Pembayaran : Dapat ditransfer ke :

Rekening BCA A/C: 215.021.3188 atas nama PT. Endo Indonesia

Rekening Mandiri A/C: 142.00.3636.777.8 atas nama PT. Endo Indonesia

(Bukti transfer dan pembayaran harap email clara.tyas@endo.id)

Demikian penawaran yang dapat kami sampaikan, besar harapan kami untuk dapat bekerjasama dan memenuhi kebutuhan rumah sakit. Apabila diperlukan informasi lebih lanjut, dapat menghubungi marketing kami di bawah ini. Terima kasih atas perhatiannya.

Hormat Kami,

PT ENDO INDONESIA



Clara Tyas Ratri

0811-3057-3636

ENDO Video Endoscopy System V3

NO.	DESCRIPTION	QTY
1	EI.EC.V3 PROCESSOR; ENDOSCOPY SYSTEM WITH SUPERIOR HIGH DEFINITION IMAGING, SMART WORK FLOW & EASY OPERATION	1
	With the HD sensor, allowing for more details to be shown in the image to improve image clarity	
	Video Output; SDI, DVI, VGA, CVBS, Y/C	
	Resolution 1080p	
	VIST1, VIST2(Variable Intelligent Staining Technology). Combined with optical and digital image processing, VIST provides perfect contrast of the image without darkening the image	
	Digital Zoom 1.4x, 1.6x, 1.8x	
	Built-in workstation offers easy access and intelligent information management. Users can manage patient data, capture images, record videos, make and print reports.	
	Patient Data >10,000 records	
	Air Feeding; 4 levels (off/L/M/H)	
	Color adjustment; Red. Blue, Saturation,	
	Structure and Edge Enhancement	
	Data Management; 500 Built-in. 5 USB, DICOM	
	AG; Yes	
	White Balance; Auto-saved	
	Contrast/Denoise; L/N/H, auto-denoise	
2	ENDO EI.LS9	1
	Main Lamp; FOUR LED lamps of different colours	
	White Light mode; High Collor rendering index and high luminous flix	
	SFI/VIST; High brightness and spesific spectra in the 2 special light modes	
	Brightness; Auto/Manual, 1-9 levels	
	Trans Illumination Mode 6-8second	
	Service time; 10000 hours (continous use)	
	Hot Plugging; Yes	
	Dimension; 473x370x125	
3	GASTROSCOPE; DISTAL END DIAMETER 9,3mm; INSTRUMENT CHANNEL DIAMETER 2,8mm; WORKING LENGTH 1050mm; TOTAL LENGTH 1390mm; FIELD OF VIEW 140 DEGREE	1
	Providing exceptional visualization of the mucosa	
	Slimmer and lighter control body with ergonomic design	

	A wide field of view angle	
	Water jet function	
	Distal end diameter (mm) : 9,3	
	Instrument channel diameter (mm) : 2,8	
	Working length (mm) : 1050	
	Total length (mm) : 1390	
	Field of view (degree) : 140	
	Observation range (mm) : 3-100	
	Bending angle : UP 210°/ DOWN 90°, LEFT 100°/ RIGHT 100°	
	Outsanding Manoeuvrability	
4	COLONOSCOPE; DISTAL END DIAMETER 12mm; INSTRUMENT CHANNEL DIAMETER 3,8mm; WORKING LENGTH 1680mm; TOTAL LENGTH 1980mm; FIELD OF VIEW 140 DEGREE	1
	High resolution and clearer image quality	
	Extra wide 140 degree field of view	
	The ergonomically designed control body has greatly increased the operability	
	Water jet function	
	Distal end diameter (mm) : 12	
	Forceps channel diameter (mm) : 3,8	
	Working length (mm) : 1680	
	Total length (mm) : 1980	
	Field of view (degree) : 140	
	Observation range (mm) : 3-100	
	Bending angle : UP 180°/ DOWN 180°, LEFT 160°/ RIGHT 160°	
	Outsanding Manoeuvrability	
5	ACCESSORIES & COMPLEMENTARY	
	MONITOR MEDICAL GRADE 26 Full HD Resolution	1
	SMART TWO" KEYBOARD	1
	Mouthpiece	5
	BIOPSY FORCEPS FOR GASTROSCOPE	5
	BIOPSY FORCEPS FOR COLONOSCOPE	5
	CLEANING BRUSH FOR GASTROSCOPE	1
	CLEANING BRUSH FOR COLONOSCOPE	1

CHANNEL PLUG FOR ENDOSCOPE	1
LEAKAGE TESTER	1
ENDOSCOPY TROLLEY + COVER	1
STABILISATOR	1
PRINTER	1
COMPRESSOR	1
ENDO Pump for Colonoscopy & Gastroscopy EI.PUMP	1
ENDO Suction Pump EI.SC7	1

ENDO Video Endoscope System EI.V2

NO.	DESCRIPTION	QTY
1	EI.EC.V2 VIDEO PROCESSOR; ENDOSCOPY SYSTEM WITH SUPERIOR HIGH DEFINITION IMAGING, VIST (VARIABLE INTELLIGENT STAINING TECHNOLOGY) AND SMART WORK FLOW & EASY OPERATION HD-500 image processor provides excellent contrast and high definition image for endoscopic examinations With the HD sensor, allowing for more details to be shown in the image to improve image clarity With VIST, gastric and rectal structures on mucous and submucosa and the distribution of vessels can be enhanced Built-in workstation offers easy access and intelligent information management. Users can manage patient data, capture images, record videos, make and print reports. (CHb) Concentration Hemoglobin Imaging Structure and Edge Enhancement Four customer-defined buttons on the control section of the endoscope handle can be preset as any main function Smart two in One keyboard; allows for doctors to use smart short-cut functional buttons and to easily manage patient information and edit reports Talented image dan video management solutions Large storage and Smart Backup Flexible reporting system Intelligent built-in workstation	1
2	EI.LS.V2; OVER 5,900K COLOR TEMPERATURE GIVE A BRILLIANT COLOR REPRODUCTION Equipped with an intelligent auto/manual lighgt control mode Main Lamp : White-LED 70W; Violet LED 50W VIST (Versalite Intelligent Staining Technology). Combined with optical and digital image processing, VIST provides perfect contrast of the image without darkening the image Average working hours : ≥10.000 hours Color Temperature : 5000K-7000K CRI : ≥90 Light control : Auto/ Manual Lightness control Air Pump : Off, Low, Mid, High Transparency : On, Off Optical Filter : On, Off Power consumption : 160 VA Frequency : 50Hz/ 60 Hz	1
3	GASTROSCOPE; DISTAL END DIAMETER 9,3mm; INSTRUMENT CHANNEL DIAMETER 2,8mm; WORKING LENGTH 1050mm; TOTAL LENGTH 1390mm; FIELD OF VIEW 140 DEGREE Providing exceptional visualization of the mucosa Slimmer and lighter control body with ergonomic design A wide field of view angle	1

	Water jet function	
	Distal end diameter (mm) : 9,3	
	Instrument channel diameter (mm) : 2,8	
	Working length (mm) : 1050	
	Total length (mm) : 1390	
	Field of view (degree) : 140	
	Observation range (mm) : 3-100	
	Bending angle : UP 210°/ DOWN 90°, LEFT 100°/ RIGHT 100°	
4	COLONOSCOPE; DISTAL END DIAMETER 12mm; INSTRUMENT CHANNEL DIAMETER 3,8mm; WORKING LENGTH 1680mm; TOTAL LENGTH 1980mm; FIELD OF VIEW 140 DEGREE	1
	High resolution and clearer image quality	
	Extra wide 140 degree field of view	
	The ergonomically designed control body has greatly increased the operability	
	Water jet function	
	Distal end diameter (mm) : 12	
	Forceps channel diameter (mm) : 3,8	
	Working length (mm) : 1680	
	Total length (mm) : 1980	
	Field of view (degree) : 140	
	Observation range (mm) : 3-100	
	Bending angle : UP 180°/ DOWN 180°, LEFT 160°/ RIGHT 160°	
5	ACCESSORIES & COMPLEMENTARY	
	MONITOR MEDICAL GRADE 24"	1
	SMART TWO" KEYBOARD	1
	Mouthpiece	5
	BIOPSY FORCEPS FOR GASTROSCOPE	5
	BIOPSY FORCEPS FOR COLONOSCOPE	5
	CLEANING BRUSH FOR GASTROSCOPE	1
	CLEANING BRUSH FOR COLONOSCOPE	1
	CHANNEL PLUG FOR ENDOSCOPE	1
	LEAKAGE TESTER	1
	ENDOSCOPY TROLLEY + COVER	1
	STABILISATOR	1
	PRINTER	1
	COMPRESSOR	1
	ENDO Pump for Colonoscopy & Gastroscopy EI.PUMP	1
	ENDO Suction Pump EI.SC7	1

ENDO Video Endoscope System EI.V7

NO.	DESCRIPTION	QTY
1	EI.EC.V7 VIDEO PROCESSOR; ENDOSCOPY SYSTEM WITH SUPERIOR HIGH DEFINITION IMAGING, SMART WORK FLOW & EASY OPERATION EI.EC.V7 image processor provides excellent contrast and high definition image for endoscopic examinations With the HD sensor, allowing for more details to be shown in the image to improve image clarity Built-in workstation offers easy access and intelligent information management. Users can manage patient data, capture images, record videos, make and print reports. (CHb) Concentration Hemoglobin Imaging Structure and Edge Enhancement Four customer-defined buttons on the control section of the endoscope handle can be preset as any main function Smart two in One keyboard; allows for doctors to use smart short-cut functional buttons and to easily manage patient information and edit reports Talented image dan video management solutions Large storage and Smart Backup 500GB Intelligent medical record management Endoscope identification DVI port; VGA port; S-Video port; CVBS port; USB port Hot plugging of endoscope Dimensions 455 mm (W) x 124 mm (H) x 370 mm (D), 9.5Kg	1
2	EI.LS.V7; OVER 5,900K COLOR TEMPERATURE GIVE A BRILLIANT COLOR REPRODUCTION Equipped with an intelligent auto/manual lighgt control mode Main Lamp : 50W LED Spare Lamp : 14W LED Average working hours : ≥50.000 hours Color Temperature : 5000K-7000K CRI : ≥90 Light control : Auto/ Manual Lightness control Air Pump : Off, Low, Mid, High Transparency : On, Off Pressure : 45kPa - 65kPa Voltage : 100V~240V AC Frequency : 50Hz/ 60 Hz	1

3	GASTROSCOPE; DISTAL END DIAMETER 9,3mm; INSTRUMENT CHANNEL DIAMETER 2,8mm; WORKING LENGTH 1050mm; TOTAL LENGTH 1390mm; FIELD OF VIEW 140 DEGREE	1
	Providing exceptional visualization of the mucosa	
	Slimmer and lighter control body with ergonomic design	
	A wide field of view angle	
	Water jet function	
	Distal end diameter (mm) : 9,3	
	Instrument channel diameter (mm) : 2,8	
	Working length (mm) : 1050	
	Total length (mm) : 1390	
	Field of view (degree) : 140	
	Observation range (mm) : 3-100	
	Bending angle : UP 210°/ DOWN 90°, LEFT 100°/ RIGHT 100°	
4	COLONOSCOPE; DISTAL END DIAMETER 12mm; INSTRUMENT CHANNEL DIAMETER 3,8mm; WORKING LENGTH 1680mm; TOTAL LENGTH 1980mm; FIELD OF VIEW 140 DEGREE	1
	High resolution and clearer image quality	
	Extra wide 140 degree field of view	
	The ergonomically designed control body has greatly increased the operability	
	Water jet function	
	Distal end diameter (mm) : 12	
	Forceps channel diameter (mm) : 3,8	
	Working length (mm) : 1680	
	Total length (mm) : 1980	
	Field of view (degree) : 140	
	Observation range (mm) : 3-100	
	Bending angle : UP 180°/ DOWN 180°, LEFT 160°/ RIGHT 160°	
5	ACCESSORIES & COMPLEMENTARY	
	MONITOR MEDICAL GRADE 21.5" Full HD Resolution	
	SMART TWO" KEYBOARD	
	Mouthpiece	
	BIOPSY FORCEPS FOR GASTROSCOPE	
	BIOPSY FORCEPS FOR COLONOSCOPE	

CLEANING BRUSH FOR GASTROSCOPE	1
CLEANING BRUSH FOR COLONOSCOPE	1
CHANNEL PLUG FOR ENDOSCOPE	1
LEAKAGE TESTER	1
ENDOSCOPY TROLLEY + COVER	1
STABILISATOR	1
PRINTER	1
COMPRESSOR	1
ENDO Pump for Colonoscopy & Gastroscopy EI.PUMP	1
ENDO Suction Pump EI.SC7	1

ENDO



ENDO EI.EC.v2 & EI.LS.v2

Video Endoscope System & Light Source


<https://endo.id>

The timeless one stop, cost-effective solution is evolving with your needs. Today, you are able to get the best out of the ENDO EI.EC.V2 image processor as we upgrade the system's configuration. Under the lasting luminance from SENDO EI.LS.V2 2-LED Light Source, you will be pleasantly surprised by your familiar video endoscopy system.




Leverage all light modes

VIST¹ is a chromoendoscopy technology combining optical and digital image processing. At the touch of a button, it provides a clear and bright image contrast, that allows a better vision of mucosa and blood vessels in the GI tract. VIST makes it easier to detect and characterize the lesions.



White light image




VIST image

Under both White Light and VIST mode, ENDO EI.EC.V2 image processor yields even more details with the 2-LED light source, and no longer sacrifices the brightness in the far field. The image quality is greatly improved.²

Whether it is a routine procedure or a complex operation, there is a benefit of having the 500 Series scopes. Their CMOS image sensor enables multi-stage enhancement in color, edge and structure, resulting in an optimal display through an internal processing.

Save more cost

The life span of the 2-LED light source is above 10,000 hours, resulting in reduced maintenance costs.




Brilliant workflow

The ENSO EI.EC.V2 features a built-in workstation with 500GB of hard disk. Storing and archiving what you have captured in a single step. It lets you access patient data management systems, generate, edit and print personalized reports smoothly, without the need of any external software.



Ergonomic scope handle design

Comfortable to hold, the 500 series scope handle adopts an ergonomic design and an intuitive layout. It helps to reduce the fatigue in longer interventions.



Outstanding manoeuvrability

Short bending radius and gradual stiffness allows to intubate, reach and observe efficiently. The large instrument channel facilitates the therapeutic procedures.



Equipped with water jet

All the scopes are equipped with the water jet function. It washes off any mucus or blood, and assures a clear view at all times.



Optimized insertion experience

The gradual stiffness makes the insertion experience easier and smoother, at any section of the GI luminal tract.

ENDO EI.EC.V2 IMAGE PROCESSOR

Station	Separated
Video output	SDI, DVI, VGA, CVBS, S-video
Resolution	960P
Digital Zoom	1.4x, 1.6x, 1.8x
Patient Data	≥ 10,000 records
IRIS	Peak/Avg/Auto
Color adjustment	Red, Blue, Chroma, ±15 levels
Air feeding	4 levels (Off/L/M/H)
AGC	Yes
Data Management	500G Built-in, 4 USB, DICOM
White Balance	Auto-saved
Image enhancement	Structure enhancement, Edge enhancement
Contrast/denoise	L/M/H, auto-denoise
Dimension (mm)	480x370x125

ENDO EI.LS.V2 LIGHT SOURCE

Main Lamp	W-LED, 70W, violet LED, 50W
White light mode	High color rendering index and high luminous flux
VIST Chromoendoscopy	Versatile Intelligent Staining Technology
Brightness	Auto/Manual, 1-19 levels
Transillumination mode	6-8 seconds
Service time	10000 hours (Continuous use)
Hot plugging	Yes
Dimension (mm)	473x370x154

SERIES SCOPE

	Gastrosopes		Colonoscopes			
	EG-500	EG-500L	EC-500	EC-500T	EC-500L	EC-500L/T
Field of view	140°	140°	140°	140°	140°	140°
Depth of view	3-100mm	3-100mm	3-100mm	3-100mm	3-100mm	3-100mm
Distal End Dia.	9.3mm	9.8mm	12mm	12mm	12.9mm	12.9mm
Insertion Dia.	9.3mm	9.8mm	12.5mm	12.5mm	12.9mm	12.9mm
Instrument Dia.	2.8mm	3.2mm	3.8mm	3.8mm	4.2mm	4.2mm
Angulation	U210°, D90°, L/R100°		U/D180°, L/R160°			
Working length	1050mm	1050mm	1350mm	1700mm	1350mm	1700mm
Water Jet	Yes	Yes	Yes	Yes	Yes	Yes

ENDO

 Jalan Raya Menganti 14 Kedurus, Surabaya 60223
Jawa Timur, Indonesia

 62-31-7673636

 62-31-7673737

 0800 . 177 . ENDO (3636)

 info@endo.id

 <https://endo.id>

ENDO



ENDO EI.EC.v7

EI.LS.v7

Video Endoscope System & Light Source

<https://endo.id>

Innovation rewrites the standard. With that in mind, ENDO Indonesia presents you the ENDO EI.EC.v7, a brand new high definition endoscopy system designed specifically for screening and diagnosis with an excellent budget friendliness.





CMOS

With the CMOS image sensor, ENDO EI.EC.v7 gives high quality images with low noise levels, which leads to improved clarity and optimized vision in not only routine procedures but also in complex interventions.



Smart Workflow

To increase efficiency, patient data management is accessible with the press of button. The smart workflow of the built-in workstation of ENDO EI.EC.v7 helps doctors to edit reports, record videos and manage patient files with ease.



Clean and Clear

The waterjet function of the scopes helps to maintain a clear endoscopic view by washing away blood or mucous.



ENDO EI.EC.v7 Video Endoscope System

HD video outputs	DVI,VGA
Standard video outputs	S-Video, CVBS
HD video recording	Built-in workstation, timely record and review
Zoom / PBP	Yes / Yes
Dimension (L x W x H / weight)	455mm x 370mm x 124mm / 9.5 kg

ENDO EI.LS.v7 Light Source

Main lamp	50W LED
Transparency function	7 seconds maximum brightness flicker
Lifetime indicator	Yes
Hot plugging	Endoscope switch without powering off
Dimension (L x W x H / weight)	460mm x 370mm x 154mm / 11kg

430 Series Scopes

Parameter	Gastroscope		Colonoscope			
	EG-430	EG-430L	EC-430	EC-430T	EC-430L	EC-430L/T
Depth of view (mm)	3-100		3-100			
Field of view	140°		140°			
Distal end diameter (mm)	9.3	9.8	12	12	12.9	12.9
Bending section diameter(mm)	9.3	9.8	12.5	12.5	12.9	12.9
Instrument channel diameter (mm)	2.8	3.2	3.8		4.2	
Angulations	UP 210°/DOWN 90° LEFT 100°/RIGHT 100°		UP 180°/DOWN 180°		LEFT 160°/RIGHT 160°	
Working length (mm)	1050		1350	1700	1350	1700
Total length (mm)	1350		1700	1950	1700	1950

ENDO

 Jalan Raya Menganti 14 Kedurus, Surabaya 60223
Jawa Timur, Indonesia

 62-31-7673636

 62-31-7673737

 0800 . 177 . ENDO (3636)

 info@endo.id

 <https://endo.id>